Kit Parcan

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Kit for canal irrigation

Properties

Although a careful mechanical canal reaming is essential, a chemical treatment is necessary for completion.

One of the procedures of this treatment consists in inserting into the canal a **relatively important amount of alkaline reaction antiseptic viz, sodium hypochlorite** which, in addition to its bactericidal potency, has the capability to desintegrate necrosed tissues and to soak off the debris stuck on the canal walls.

- This irrigation should not be irritant for the apical area. That is the reason why the sodium hypochlorite solution recommended is diluted.
- So as to be efficient, the amount of liquid to be used should be relatively important (10 to 20 ml) in order that all the debris might be soaked off and dragged away.
- So as to avoid injecting hypochlorite solution into periapical tissues, it is recommended to have **an appropriate needle** in order that the liquid injected may not, in any case, penetrate into the periapex.

The proposed needle has two apertures at its tip : one orifice at the end and one lateral perforation. This lateral perforation presents the following advantage : in the case where the irrigation needle tip is blocked in a narrow part of the canal, the liquid could be injected beyond the apex because of the pressure exerted on it. The lateral perforation allows this liquid to escape and reach wider areas of the canal and make its way towards the pulp chamber with efficient activation.

Instructions for use

Pump the required amount of hypochlorite from the flask. Then adjust the needle onto the syringe and carry out the canal irrigation.

As with any other endodontic instruments, the needle should be carefully inserted into the canal in order to carry out the treatment without injuring the periapex.

It must be calculated that the point should remain at about 1.5 mm from the apex, considering that the irrigation liquid has a low surface tension and a low pressure. Before inserting the needle into the canal, it can be bent according to a desired length. A perforated disk can also be slipped onto the needle so as to gauge the length of the needle that should penetrate into the canal.

Warnings

KIT PARCAN instruments must not be sterilized by heat.

Presentation

The kit contains :

- one 10 ml syringe
- 2 perforators (yellow tips)
- 10 blunt pointed needles, with an orifice at the end and a lateral perforation (orange tips)

For professional use only in the practice of dentistry.