

Dental Syringes

Septodont manual aspirating syringes



DESCRIPTION

High quality, Type CW winged, manual aspirating syringe, with harpoon, available in two grip sizes.

Both sizes accept standard 1.7 mL and 1.8 mL dental cartridges.

PRESENTATION

Each, two sizes available, large or small grip.
Please specify when ordering

Item	Description	Packaging
11049	Septodont aspirating syringe large (gold)	Each
11050	Septodont aspirating syringe small (blue)	Each



Septodont of Canada, Inc.

Distributors of Septodont and Novocol Products
Cambridge, ON (800) 647-0643
www.septodont.ca